
## 样品信息

: 2024/09/30 15:53:52 : 18

: D:\DATA\白酒-手性\全扫定性\2024\0929\1\18.qgd : D:\DATA\白酒-手性\全扫定性\2024\0929\18.qgd : D:\DATA\白酒-手性\method\酒-Scan-240912-乙醇.qgm : D:\DATA\白酒-手性\method\酒-Scan-240912-乙醇.qgm



			峰表 TIC						
峰号	保留时间	I.Time	F. Time	峰面积	峰面积%	峰高	CAS#	名称	
1	2. 056	2.000	2.240	656870	0.0209	96695	75-07-0	乙醛	
$\bar{2}$	2. 582	2. 530	2.695	291044	0.0093	68605	67-64-1	丙酮	
3	3. 234	3. 170	3. 445	4807125	0. 1531	1197578	141-78-6	乙酸乙酯	
4	4. 713	3. 865	4. 930	1329631195	42. 3427	33518611	64-17-5	乙醇	
5	5. 023	4. 950	5. 100	229615	0.0073	53440	107-87-9	2-戊酮	
6	6. 314	6. 260	6.375	67332	0.0021	21378	78-92-2	仲丁醇	
7	6. 457	6.395	6.545	739067	0.0235	212575	105-54-4	丁酸乙酯	
8	6. 777	6.650	6.870	621044	0.0198	128229	71-23-8	正丙醇	
9	10. 198	6.910	10.630	1780776305	56.7097	12753131	7732-18-5	水	
10	11.077	10.990	11.170	975275	0.0311	304633	71-36-3	正丁醇	
11	13. 304	13. 130	13.480	499585	0.0159	114190	123-51-3	异戊醇	
12	13.632	13.480	13.730	10076240	0.3209	2393010	123-66-0	己酸乙酯	
13	15.008	14.945	15.060	44824	0.0014	13816	71 - 41 - 0	正戊醇	
14	16. 507	16.420	16.540	23917	0.0008	6648	513-86-0	3-羟基-2-丁酮	
15	18.007	17.915	18.095	68513	0.0022	21145	106-30-9	庚酸乙酯	
16	18.928	18.780	19. 105	3440777	0.1096	843168	687-47-8	L-乳酸乙酯	
17	19. 280	19. 195	19.380	301457	0.0096	83859	111-27-3	正己醇	
18	21.437	21.360	21.530	42323	0.0013	10825	626-82-4	己酸丁酯	
19	22.399	22.305	22.485	71625	0.0023	19484	106-32-1	辛酸乙酯	
20	23.808	23.680	24.055	2742720	0.0873	511447	64 - 19 - 7	乙酸	
21	27.076	26.985	27.120	28373	0.0009	8639	25601-36-9	2-羟基乙基己酸酯	
22	27.165	27.120	27.235	44990	0.0014	10481	513-85-9	2,3-丁二醇	
23	27. 291	27.235	27.385	36342	0.0012	7198	79-09-4	丙酸	
24	28. 414	28.345	28.500	31649	0.0010	8664	79-31-2	异丁酸	
25	29.272	29. 180	29.370	32964	0.0010	6911	57-55-6	丙二醇	
26	30. 726	30.645	30.865	585617	0.0186	143195	107-92-6	丙二醇甲醚	
27	32.252	32.195	32.360	42078	0.0013	9573	503-74-2	异戊酸	
28	34. 788	34.695	34.925	106348	0.0034	25352	109 - 52 - 4	正戊酸	
29	38. 593	38. 480	38.805	3042554	0.0969	659622	142-62-1	正己酸	
30	40.552	40.475	40.625	35300	0.0011	8030	60-12-8	苯乙醇	
31	42. 230	42.170	42.320	20604	0.0007	5529	111-14-8	庚酸 辛酸	
32	44. 364	44.305	44.430	35406	0.0011	15662	124-07-2	辛酸	
33	45.874	45.815	45.915	14311	0.0005	5423	4511-42-6	L-乳酸	
				3140163389	100.0000	53286746			